



Fact Sheet (5)

freeflex®: The Innovative IV-Bag

As leader in infusion and nutrition therapy, Fresenius Kabi is one of the largest and most capable manufacturers of IV-bags in the world. Our expertise in this area is continuously challenged by our drive to develop even better products than already available: safer for the patient, more convenient to handle at the hospitals, more environmentally friendly.

Freeflex®, Fresenius Kabi's innovation bag, represents an improvement in all of these areas! Freeflex® is an IV-bag of unparalleled qualities using state-of-the-art, fully integrated manufacturing processes. Freeflex® is an entirely PVC-free bag-system that shows out-standing drug compatibility at an extended pH interval and is extremely user-friendly.

The development of Freeflex® puts Fresenius Kabi Product Partnering in an even better position to tailor our service exactly to the needs and wishes of our clients!

freeflex®: The advantages

- Easy Handling
- Enhanced Safety
- Excellent Drug Compatibility
- Environmentally Friendly

freeflex®: In a nutshell

Container sizes ranging from 50ml up to 1000ml

Strong and resistant central hanger

High strength of welded seam enabling the use of pressure cuff

Safe and simple reconstitution avoiding any airborne contamination using the Freeflex® transfer device

Multi-layer PP primary bag film, free of PVC, plasticizers, additives or glues, with excellent (!) drug compatibility at a large pH range, low gas and water (!) vapour permeability, high transparency, and resistant to sterilization at 121 °C

Easy to remove, peelable overwrap

Tamper evident break-off ports with sterile chamber. Easy to identify and easy to open. Due to a special design of the ports, leakage is prevented, even after the removal of the administration set.

freeflex®: Award winning

Freeflex Wins Top French Award for Innovation

Fresenius Kabi's revolutionary non-PVC intravenous fluid and drug combination delivery system, Freeflex®, now available in the UK, has been awarded the top accolade for innovation by leading independent hospital physicians and pharmacists in France.

Taken from www.inside-hospitals.co.uk